


MgFeO XRD-2D combined analysis using MAUD

macro-strain
texture
crystallite size
micro-strain


a tutorial transcription from the XRD school in Trento December 2018




MAUD image 2D calibration

1) Instrument calibration:

Dataset : DataFileSet_x: Edit



Instrument Edit



Instrument Name

Diffraction Instrument

File Edit

Instrument name: 2DImage_ALS25keV

Incident intensity: 2.0

Intensity calibration: none cal Options

Angular calibration: Flat Image Transmission Options

Geometry: Image 2D Options

Measurement: 2Theta Options

Source: Synchrotron Options

Detector: Scintillation Options

Instrument Broadening: Caglioti PV Options

OK

Angular calibration: Flat Image Transmission

Detector distance: 285 mm

Flat Image in transmission/reflection angular calibration

File Edit

Obsolete, use Inclined Reflection Image with 2theta=0 (transmission) and 2theta=180 (reflection)

Calibration parameters Integration setting

Detector distance: 285 Center x error: 0

Center y error: 0 Tilting error x: 0

Tilting error y: 0 Ratio width/height pixels: 1.0

Cancel OK

Geometry: Image 2D

Measurement: 2Theta (means no scan, as opposed to Theta-2Theta)

Source: Synchrotron: wavelength: 0.49594 Å (25 keV)

Radiation

File Edit

+ add term

- remove term

Added_by_d...

Wavelength: 0.49594


Weight: 1.0

Cancel OK


Detector: Scintillation (only change for XRF)

Instrument Broadening: Cagliotti PV


Asymmetry: 0, 0 (Fixed); Use reciprocal value: check




HWHM: 0.0025, 0, 0 (Fixed)



2) DataFiles: From Images: opens ImageJ 1.51r26 (last version working with MAUD)




File Open: LaB6_25keV_003.tiff16



Image>Adjust>Brightness/Contrast:


Set: Min=0, Max=255



Check if the image has size in pixels and mm: we need the pixel size
Image>Properties:

Create a mask to avoid all pixel with value = 0 or value < 0


Use “Oval” selection tool



and adjust a circle a little bit inside the exposed area:

then,

Plugins>Maud Plugins>Multi-Spectra from normal transmission/Reflection Image



Distance: 285 mm

center X: 182.3 mm

center Y: 162.7 mm

Tracker circle radius: 17.0 deg (2tht for LaB6)

[Update]

Integrate from 0 deg to 360 deg, 18 spectra


For capillary: omega=0, chi=0, phi=0

(for slabs: explanation in 20181129_143839-46-cont.mp4 at t=00:28:30)

[OK] , wait for integration and save with “.esg” extension


Close ImageJ, Don't save image changes

and get in MAUD





then check:

General>Instrument>Edit>Angular Calibration>Flat Image Transmission>Options>
>Integration setting: center X, center Y



close, and you get:





and in [Plot 2D] one can see the beam-stopper in spectrum #9

DataFiles> file .esg(9) > Enable : check disable

Also observe small deviation due to centering

Datasets> General > Range of Computation> Min in data units: 6 and !!! Max : 40
(if you leave Max: 0, the program will stuck !!!)

Phases>Load from CIF: LaB6.cif Lattice value: 4.15689 A

Fix lattice value and refine detector distance and centering (also tilt), and instrument broadening

To fit instrument broadening:

Phase>Edit>Microstructure> Line broadening model>None


or

Line broadening model>Delft

and

Size-strain model>Isotropic>Options> 0, 0

Then, [Compute spectra]



We observe that the distance is not good: change to 350 mm and calculate again:

and start to refine:

Background: normal one ~ 100

Intensity: (instrument edit) : decrease 20-40 times

Analysis> Parameter list> Expand All:

and [Free Backgrounds]

instrument> incident intensity: refine! (`_pd_proc_intensity_incident`)


Flat Image Transmission: refine!:

distance, center x,y, tilt x,y

Image 2D: polarisation effects for highly polarized beam

(in the vertical direction, or horizontal at high angles; otherwise: 0,0)

[Fit]



Try to fit the grain-iness using arbitrary texture (not actually texture but arbitrary intensity)

Phase> Edit>Advanced model>Texture>Arbitrary tex

before Fit, fix incident intensity !!!!

[Fit] again

Next, zoom a peak and observe the asymmetry:

Instrument>Cagliotti>Asymmetry 1st parameter refine

and Truncation angle: set negative! -0.4

means, asymmetry increase toward higher angles (for 2D detectors)

[Fit]

Refine 2nd asymmetry parameter

Refine HWHM Cagliotti parameters

[Fit] :

Refine Gaussian Cagliotti parameters

[Fit] :

Dataset>Edit>Exclude region>Add> from 15.5 to 16.0

Analysis>Parameter List>Fix All Parameters

Dataset>Instrument>Save>LaB6_ALS_image2D.ins

2) MgFeO Magnesiumwuestite high pressure

New Analysis

Dataset>Instrument>Import>LaB6_ALS_image2D.ins

Dataset>DataFiles>From Image>ImageJ

ImageJ: File>Open>MgFeO_25keV_006.tiff16

ImageJ: Image>Adjust>Brightness/Contrast: [Set] min=0, max=255


ImageJ: Toolbar: Oval: Select exposed area


ImageJ: Plugins>Maud Plugins>Multi-Spectra from normal transmission/Reflection Image

Number of spectra: 72 (every 5 deg)

[ok]

save *.esg file






Back in MAUD

Zoom in at first line: change dataset range: min=11.5, max=50


Phase>Import>MgFeO.cif

[Calculate]




Shift is caused by pressure: 40 Gpa

Adjust lattice parameter: from 4.24 to ~ 4.0 Å




Also in Plot 1D



Start analysing:

Background, and scale parameters



Remove beam-stopper: `esg(35) disable`

NEXT: Stress model

Phases>Edit>Advanced Models:

Strain: Moment Pole Stress > Options

> stiffness matrix at 40 gPa for MgFeO

C11: 578.25

C12: 161.9

C44: 141.4

> hydrostatic stress (Macrostress)


S11=S22

S33= -2*S11

neglect shear: S12, S13, S23

“One question ... From where do we get these value?”

“From literature ...”




Differential stress:

$$\sigma_{ij} = \begin{bmatrix} \sigma_p & 0 & 0 \\ 0 & \sigma_p & 0 \\ 0 & 0 & \sigma_p \end{bmatrix} + \begin{bmatrix} -t/3 & 0 & 0 \\ 0 & -t/3 & 0 \\ 0 & 0 & 2t/3 \end{bmatrix} = \sigma_p + D_{ij}$$

Let's say pressure along S22 = 3Gpa (same units as stiffness matrix) and Right-click: "refine"
Then S11 is minus half of S22


right click S11> select "Equal to"



[Set bound]

Do same for S33

[OK] and open from Options again, and check:




Stress/strain model: BulkPathGEO (geometrical mean path)

Refine also:

- > lattice parameters
- > crystallite size and micro-strain

Before, just check: >Compute >Nothing !?

Verify: >Graphic>Texture plot>Pole Figure Coverage> (2,0,0) and (2,2,2) > Plot



It's OK: the reason is the stress is along S11 and not S33 . Let's change:

S11 = -3 Gpa (Refine)

S22 (set equal to) S11 : * (-0.5) + (0)

S33 (set equal to) S11 : * (-0.5) + (0)

Moment pole figures options panel

File Edit

Macrostress11 : -3 Macrostress22 : 1.5 Stress/strain model: BulkPathGEO

Macrostress33 : 1.5 Macrostress23 : 0 Weight (Voigt-Reuss): 0.5


Macrostress13 : 0 Macrostress12 : 0 ☐ Use texture ODF

Stiffness matrix

578.25	161.9	161.9	0	0	0
-	578.25	161.9	0	0	0
-	-	578.25	0	0	0
-	-	-	141.4	0	0
-	-	-	-	141.4	0
-	-	-	-	-	141.4

Cancel OK


Calculate: OK



OK, let's refine:

- > Background
- > Incident intensity
- ! Do NOT refine detector tilt
- ! 1st, keep detector centering fixed
- > crystallite size and microstrain
- [Fit]

$$\begin{array}{c}
 | \\
 | \\
 \text{vS11} \\
 / \text{-----} \rightarrow \text{S33} \\
 / \\
 | \text{S22}
 \end{array}$$



Now we need Texture

Phase>Edit>Advanced Models>**Texture**>Standard Functions>Options
 Fiber or Spherical


Here: **Fiber along S11**

add component

ThetaY = 90

PhiY = 0


PhiH = c axis




ThetaH: the crystallographic plane along the fiber component: 100 = 001

For 001: ThetaH=0, PhiH=0


FWHM: angular spread of the fiber: 30 deg quite sharp, but we will refine



Check if direction is OK:
Calculate




Graphic>Texture Plot>Select : 111, 200, 220
Reconstructed intensity: check > Plot!




See fiber along 200 = 002

Refine Fiber texture: FWHM and gaussian parameter: Now OK



If we do not know about the texture fiber, and we start with with a model without imposing any symmetry.


Phase>Edit>Advanced Model> Texture>E-WIMV > Options




symmetry: none

ODF resolution: 15 deg (start rough and later refine mor fine ~ 7.5 deg)
 OK and refine!


Graphic>Texture Plot>Select : 111, 200, 220
 Reconstructed intensity: check > Plot!



The fiber can be recognized; let's use higher resolution: 5 deg
 Refine!



And inverse pole figures:



111

100

110